

Program for evaluation of the 3ω results

August 11, 2016

Contents

| | | |
|----------|--|----------|
| 1 | General information | 2 |
| 1.1 | Names of the files | 2 |
| 1.2 | Structure of the files | 2 |
| 2 | Subprogram HeaterWidth | 4 |
| 2.1 | Single samples: | 4 |
| 2.2 | Series of samples: | 4 |
| 3 | Subprogram Dialog3ω | 4 |
| 4 | Governing equations | 5 |

To evaluate the 3ω data run the file *Dialog3w.m*.

1 General information

1.1 Names of the files

- file with 3ω results: SampleName_PowermW_U3w.dat (e.g. DSO-NSTO-500p_15mW_U3w.dat)
- file with dRdT data: SampleName_dRdT.dat (e.g. DSO-NSTO-500p_dRdT.dat)
- file with current, voltage and relative width of the heater: SampleName_PowermW_VIP.dat (e.g. DSO-NSTO-500p_25mW_VIP.dat)
- images: SampleName_0x.png (e.g. DSO-NSTO-500p_01.png)

1.2 Structure of the files

1. File with 3ω results (data in μV):

| f (Hz) | Sample1 | Sample2 | Sample3 | Sample4 | // structure |
|--------|-------------|---------|---------|---------------|--------------|
| | DSO-NSTO | empty | empty | DSO-Substrate | |
| 100 | 1211.77002 | NaN | NaN | 1133.449951 | |
| 116 | 1164.210083 | NaN | NaN | 1090.650024 | |
| 135 | 1118.910034 | NaN | NaN | 1047.97998 | |
| 157 | 1073.01001 | NaN | NaN | 1006.369934 | |
| 183 | 1027.710083 | NaN | NaN | 963.4530029 | |
| 212 | 984.5569458 | NaN | NaN | 923.756958 | |
| 247 | 940.5650024 | NaN | NaN | 882.1519775 | |
| 287 | 897.6500244 | NaN | NaN | 842.4559937 | |
| 333 | 855.9260254 | NaN | NaN | 803.473999 | |
| 387 | 814.7990112 | NaN | NaN | 764.7310181 | |
| 450 | 774.1490479 | NaN | NaN | 726.4650269 | |
| 523 | 734.690979 | NaN | NaN | 689.2720337 | |
| 608 | 696.0670166 | NaN | NaN | 652.9130249 | |
| 707 | 658.3969727 | NaN | NaN | 617.5079956 | |
| 822 | 621.6799927 | NaN | NaN | 583.2949829 | |
| 955 | 586.7520142 | NaN | NaN | 550.1550293 | |
| 1110 | 552.7770386 | NaN | NaN | 518.2069702 | |
| 1291 | 520.1139526 | NaN | NaN | 487.4490051 | |
| 1500 | 489 | NaN | NaN | 458.1229858 | |
| 1744 | 458.9580078 | NaN | NaN | 429.9899902 | |
| 2027 | 430.3469849 | NaN | NaN | 403.0480042 | |
| 2356 | 403.2870178 | NaN | NaN | 377.5379944 | |
| 2739 | 377.5379944 | NaN | NaN | 353.4570007 | |
| 3183 | 353.3379822 | NaN | NaN | 330.8079834 | |
| 3700 | 330.4500122 | NaN | NaN | 309.4689941 | |

if you have 2 samples instead of 4, it is necessary to write NaNs in 2 columns so these columns will not be taken into account while calculation.

Important!

- Column 1: frequency (Hz)
- Columns 2-3: thick films ($U_{3\omega}(\mu V)$)
- Columns 4-5: reference films ($U_{3\omega}(\mu V)$)

2. File with dRdT data (data must be in Ohm!):

| Pt100 | Sample1 | Pt100 | Sample2 | Pt100 | Sample3 | Pt100 | Sample4 |
|--------|---------|-------|---------|-------|---------|--------|-----------|
| Pt100 | S1 | empty | empty | empty | empty | Pt100 | Substrate |
| 107.13 | 24.84 | NaN | NaN | NaN | NaN | 107.12 | 23.66 |
| 109.75 | 25.19 | NaN | NaN | NaN | NaN | 109.76 | 23.99 |
| 111.22 | 25.39 | NaN | NaN | NaN | NaN | 111.24 | 24.20 |
| 112.46 | 25.56 | NaN | NaN | NaN | NaN | 112.48 | 24.36 |

The data saved in the files from dRdT data acquisition is measured voltage in mV. The current applied during the measurement is $100 \mu A$. Therefore, to get the resistance in Ohm one has to multiply the measured data by 10. Additionally, in the SampleName_dRdT_measured.dat file the measurements for positive and negative current are saved. Therefore, it is necessary to calculate an average from both entries - first(third...) line corresponds to a measurement for +I whereas second(fourth...) line is V measured for -I.

Important!

- Columns 1-4: thick films ($R (\Omega)$)
- Columns 5-8: reference films ($R (\Omega)$)

3. File with current and voltage used during the 3ω measurement and relative width of the heater:

| Quantities | Sample1 | Sample2 | Sample3 | Sample4 | //structure |
|------------|-------------|---------|---------|-------------|-------------|
| | DSO-NSTO | empty | empty | DSO | Substrate |
| Uappl (V) | 3.971999884 | NaN | NaN | 4.427999973 | |
| Uw (V) | 0.819400012 | NaN | NaN | 0.800800025 | |
| Iw (mA) | 30.52700043 | NaN | NaN | 31.22599983 | |
| Pixel | 0.93025 | NaN | NaN | 0.94458 | |

Important!

- Column 1: names of the quantities
- Columns 2-3: thick films
- Columns 4-5: reference films

Uappl (V), Uw (V), Iw (mA) are saved during the 3ω measurement. Quantity called "pixel" is a relative heater width that is determined by the subprogram called "HeaterWidth". It is width of a heater on a given sample (in pixels) divided by the width of a heater on a mask (also in pixels). The width of mask for $20 \mu m$ structure is saved in the program.

2 Subprogram HeaterWidth

There are two main layouts of this program - designed for single samples (Layout/Single samples) and for a series of samples (Layout/series of samples). The second variant is default.

The program works in the following way: the user specifies two points on one side of the heater line that create a line. In the second step the user chooses several points (e.g. 5) on the opposite edge and the program calculates a distance between the line and these points and calculates an average. In Single samples mode this is the outcome of the analysis. In Series of samples mode the final result is an average of e.g. 5 points per image and a given number of images (e.g. 3), so in total there are 15 heater widths to be taken to the average calculation.

The outcome of the program is a relative heater width - it is width of the heater on the given sample (in pixels) divided by the width of the heater on a mask (also in pixels). The width of mask for 20 μm structure is saved in the program.

2.1 Single samples:

- Pressing on buttons Sample 1, Sample 2, etc. you choose the files you want to measure.
- Specify number of points to be measured - how many measurements do you want to take for every line and press OK button. If you choose for example 5, then a program calculates an average of the 5 measurements and gives you a mean value as an outcome.
- When you are done with choosing the points the result will show up in the small table. You can save the data pressing Save button.

2.2 Series of samples:

- Specify the folder where the data for thick and reference samples are saved.
- Filenames for images should have form SampleName_01.png. In a field next to "Sample" label you have to write SampleName and specify how many images are saved in a series. E.g. if you have S1_01.png, S1_02.png, S1_03.png and S1_04.png, you write "S1" next to label "Sample 1" and "4" next to "No. of images".
- If you have only 2 series of samples to measure that are in the same directory, you have to specify only the path to "Thick samples", however, if you have 4 series of samples, you have to specify "Ref. samples" path as well. Even if all the images are in one folder you have to specify the path twice.

3 Subprogram Dialog3 ω

- Specify the directory that contains your data. Pay attention to the filename pattern that was described at the beginning of this document (section 1.1).
- Insert your SampleName and Power - based on this information the program will be able to find your files.
- Insert the thickness of your sample - the difference in thickness between thick and the reference samples.
- Specify the dimensions of a heater structure that you use. Be careful! The program asks you for b, not 2b - if total width of your heater line is $2b = 20 \mu\text{m}$, you have to insert $b = 10 \mu\text{m}$.

- You can save your results in .png or/and .pdf files. Three files will be saved:
 - dRdT fit,
 - fit of the $U_{3\omega}$ for slope based method and the determination of substrate thermal conductivity,
 - temperature oscillations for every sample and the difference between temperature oscillations for thick and reference samples
- additionally program generates also .xlsx (Excel) file containing summary of all the calculated data.

4 Governing equations

- $P = U_{\omega} \cdot I_{\omega}$
- $R = \frac{U_{\omega}}{I_{\omega}}$
- $\frac{dR}{dT} = \frac{dR_{heater}}{dT} = \frac{dR_{Pt100}}{dT} \cdot \frac{dR_{heater}}{dR_{Pt100}} = 0.390802 \frac{\Omega}{K} \cdot \frac{dR_{heater}}{dR_{Pt100}}$
- Thermal conductivity of the substrate - calculated based on a slope method

$$\kappa_{substrate} = \frac{U_{\omega}^3}{4\pi l R^2} \frac{1}{\frac{d(U_{3\omega})}{d(\ln 2\omega)}} \frac{dR}{dT}$$

- Temperature oscillation introduced to the sample

$$\Delta T = 2 \frac{1}{\frac{dR}{dT}} R \frac{U_{3\omega}}{U_{\omega}}$$

- Considering a heat flow analogously to an electric current, one can write an analogy to the Ohm's law.

| | |
|--|---|
| $I = \frac{U}{R} = \frac{\sigma S}{L} U$ | $Q = P = \frac{\kappa_{TF} \cdot 2b \cdot l}{d_{TF}} \Delta T_{TF}$ |
| σ | κ |
| S | $2b \cdot l$ |
| L | film thickness d_{TF} |
| U | ΔT_{TF} |

where $d_{TF} = d_{thick\ film} - d_{reference}$ and $\Delta T_{TF} = \Delta T_{thick\ film} - \Delta T_{reference}$, $2b$ is the heater line width and l is its length between the voltage pads, P is the power used for a measurement. Based on the comparison presented above, thermal conductivity of the thin film may be calculated according to

$$\kappa_{TF} = \frac{P \cdot d_{TF}}{2b \cdot l \cdot \Delta T_{TF}}$$